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
# Anti-MPO (P-ANCA) ELISA

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## REFERENCES

1. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, van Es LA, van der Giessen M, van der Hem GK, The TH: Autoantibodies against neutrophils and monocytes: tool for diagnosis and marker of disease activity in Wegeners granulomatosis. *Lancet* i:425-429, 1985.
2. Gross WL, Lüdemann G, Kiefer G, Lehmann H: Anticytoplasmic antibodies in Wegeners granulomatosis. *Lancet* i:806,1986.
3. Falk RJ, Jennette JC: Immunofluorescence and ELISA determination of ANCA with description of a sub-class with anti-myeloperoxidase activity. In: Rasmussen N, Wiik A (eds): Anti neutrophil cytoplasm antibodies. Proceedings of the 1st International Workshop on ANCA (Copenhagen, Denmark, 25-26 January 1988). *APMIS* 97(suppl 6), 1989.
4. Falk RJ, Jennette JC: Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 318:1651-1657, 1988.
5. Andrassy K, Koderisch J, Waldherr R, Rufer M: Diagnostic significance of anticytoplasmic antibodies (ACPA/ ANCA) in detection of Wegeners granulomatosis and other forms of vasculitis. *Nephron* 49:257-258, 1988.
6. Specks U, Wheatley CL, McDonald TJ, Rohrbach MS, DeRemee RA: Anticytoplasmic autoantibodies in the diagnosis and follow-up of Wegeners granulomatosis. *Mayo Clin Proc* 64:28-36, 1989.
7. Nölle B, Specks U, Lüdemann J, Rohrbach MS, DeRemee RA, Gross WL: Anticytoplasmic autoantibodies: their immunodiagnostic value in Wegener granulomatosis. *Ann Intern Med* 111:28-40, 1989.
8. Lüdemann J, Csernok E, Ulmer M, Lemke H, Utecht B, Rautmann A, Gross WL: Anti-neutrophil cytoplasm antibodies in Wegener's granulomatosis: immunodiagnostic value, monoclonal antibodies and characterisation of the target antigen. In: de Leeuw PW (ed): Proceedings of the 2nd International Workshop on antineutrophil cytoplasmic antibodies (ANCA) (Noordwijkerhout, The Netherlands, 23-24 May 1989). *Neth J Med* 36:157-162, 1990.
9. Lüdemann J, Utecht B, Gross WL: Anti-neutrophil cytoplasm antibodies in Wegener's granulomatosis recognize an elastolytic enzyme. *J Exp Med* 171:357-362, 1990.
10. Kallenberg CGM: Antineutrophil cytoplasmic antibodies (ANCA) and vasculitis. *Clin Rheumatol* 9(suppl 1):132-135, 1990.
11. Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The TH, van der Hem GK, Kallenberg CGM: Prevention of relapses in Wegeners granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. *Lancet* 336:709-711, 1990.
12. Cohen Tervaert JW, Goldschmeding R, Elema JD, Limburg PC, van der Giessen M, Huitema MG, Koolen MI, Hene RJ, The TH, van der Hem GK, van dem Borne AEGK, Kallenberg CGM: Association of autoantibodies to myeloperoxidase with different forms of vasculitis. *Arthritis Rheum* 33(8):1264-1272, 1990.
13. Ulmer M, Rautmann A, Gross WL: Immunodiagnostic aspects of autoantibodies against myeloperoxidase. *Clin Nephrol* 37(4):161-168, 1992.
14. Kallenberg CGM, Mulder AHL, Cohen Tervaert JW: Antineutrophil cytoplasmic antibodies: A still-growing class of autoantibodies in inflammatory disorders. *Am J Med* 93:675-682, 1992.
15. Bini P, Gabay JE, Teitel A, Melchior M, ZHOU J, Elkon K: Antineutrophil cytoplasmic autoantibodies in Wegener's granulomatosis recognize conformational epitope(s) on proteinase 3. *J Immunol* 149(4):1409-1415, 1992.
16. Noël L-H et al.: ANCA: Diversity and clinical applications. *Advances Nephrol* 22:238-267, 1993.
17. Jennette JC et al.: Nomenclature of systemic vasculitides. Proposal of an International consensus conference. *Arthritis Rheum* 37(2):187-192, 1994.
18. Esnault VLM, Short AK, Audrain MAP, Jones SJ, Martin SJ, Skehel JM, Lockwood CM: Autoantibodies to lactoferrin and histone in systemic vasculitis identified by anti-myeloperoxidase solid phase assays. *Kidney International* 46:153-160, 1994.
19. Schnabel A, Csernok E, Isenberg DA, Mrowka C, Gross WL: Antineutrophil cytoplasmic antibodies in systemic lupus erythematosus. Prevalence, specificities, and clinical significance. *Arthritis Rheum* 38(5):633-637, 1995.
20. Forde AM, Feighery C, Jackson J: Characterisation of anti-neutrophil cytoplasmic antibody target antigens using electrophoresis and Western blotting techniques. *British Journal of Biomedical Science* 55:247-252, 1998.
21. Savige J, Gillis D, Benson E, Davies D, Esnault V, Falk RJ, Hagen EC, Jayne D, Jenette JC, Paspaliaris B, Pollock W, Pusey C, Savage COS, Silvestrini R, van der Woude F, Wieslander J, Wiik A: International consensus statement on testing and reporting of antineutrophil cytoplasmic antibodies (ANCA). *Am J Clin Pathol* 111:507-513, 1999.
22. Gross WL, Trabandt A, Reinhold-Keller E: Diagnosis and evaluation of vasculitis. *Rheumatology* 39(3):245-252, 2000.
23. Wiik A: Relevant target antigens of ANCA in primary vasculitis: which test should be used in clinical practice? In: Walport MJ (ed): Proceedings of the 9th international ANCA workshop (Groningen, The Netherlands, 12-15 April 2000). *Clin Exp Immunol* 120 (suppl 1): 9-10, 2000.
24. Boomsma MM, Stegeman CA, van der Leij MJ, Oost W, Hermans J, Kallenberg CG, Limburg PC, Cohen Tervaert JW: Prediction of relapses in Wegener's granulomatosis by measurement of antineutrophil cytoplasmic antibody levels: a prospective study. *Arthritis Rheum* 43 (9): 2025-2033, 2000.

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