




Instruction for use


Serotonin - ELISA

Enzyme Immunoassay
for the Quantitative Determination of Serotonin
in Serum, Plasma, Urine, and Cerebrospinal Fluid



REF EA602/96

 12 x 8

 2 – 8 °C

DLD Gesellschaft für Diagnostika und medizinische Geräte mbH
Adlerhorst 15 • D-22459 Hamburg • Tel +49-40-555 87 10 • Fax +49-40-555 87 111
Internet: <http://www.dld-diagnostika.de> • E-Mail: contact@dld-diagnostika.de

April 2009

Contents

| | | |
|---|------|----|
| 1. Introduction and Principle of the Test | Page | 4 |
| 2. Precautions | Page | 4 |
| 3. Storage and Stability | Page | 4 |
| 4. Contents of the Kit | Page | 5 |
| 5. Sample Collection | Page | 6 |
| 6. Preparation of Reagents and Samples | Page | 7 |
| 7. Test Performance (ELISA) | Page | 8 |
| 8. Calculation of the Results | Page | 9 |
| 9. Assay Characteristics | Page | 10 |
| 10. Literature | Page | 14 |
| Pipetting Scheme | Page | 16 |

1. Introduction and Principle of the Test

Serotonin (5-Hydroxytryptamine), a biogenic amine, is a product of the tryptophan metabolism. It is a well evaluated neurotransmitter of the central nervous system and can be found in high concentrations in the chromaffine cells of the intestinal mucosa, in the platelets and the serotonergic neurones of the brain.

Central-serotonergic neurones influence physiological functions such as sleep and the hormonal and cardio-vascular regulation. Increased serum levels can be found with malignant carcinoid, endogenous depression and schizophrenia.

The assay kit provides materials for the quantitative measurement of derivated serotonin (5-Hydroxytryptamine) in serum, plasma and urine. The derivation is performed during the preparation of the samples. By using the acylation reagent the serotonin is quantitatively derivated into N-acylserotonin.

The competitive Serotonin ELISA kit uses the microtitre plate format. Serotonin is bound to the solid phase of the microtiter plate. Acylated serotonin and solid phase bound serotonin compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. The antibody bound to the solid phase serotonin is detected by anti-rabbit/peroxidase. The substrate TMB / peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase serotonin is inversely proportional to the serotonin concentration of the sample.

2. Precautions

- For in vitro use only.
- Disposable gloves and safety glasses should be used.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy but these materials should be handled as potentially infectious.
- Some components of this kit are containing hazardous reagents. These components are marked with the adequate hazard label.

3. Storage and Stability

On arrival, store the kit at 2-8 °C. Once opened the kit is stable until its expiry date. For stability of prepared reagents refer to Preparation of Reagents.

Do not use components beyond the expiration date shown on the kit labels.

Do not mix various lots of any kit component within an individual assay.

4. Contents of the Kit

4.1 **MT-Strips** **STRIPS** 12 strips
8 wells each, break apart
precoated with serotonin

4.2 **Standards 1 - 6** **CAL 1-6** 6 vials
Each 4 ml, ready for use
Concentrations:

| Standard | 1 | 2 | 3 | 4 | 5 | 6 |
|----------|---|----|----|-----|-----|-------|
| ng/ml | 0 | 15 | 50 | 150 | 500 | 2,500 |

4.3 **Control 1 & 2** **CONTROL 1 & 2** 2 vials
Each 4 ml, ready for use
Range: see q.c. certificate

4.4 **Acylation Buffer** **ACYL-BUFF** 1 vial
3.5 ml, ready for use

4.5 **Acylation Reagent** **ACYL-REAG** 1 vial
3 ml, ready for use

4.6 **Antiserum** **AS** 1 vial
5.5 ml, ready for use, colour coded blue
Rabbit-anti-N-acylserotonin

4.7 **Enzyme Conjugate** **CONJ** 1 vial
11 ml, ready for use
Goat anti-rabbit-IgG-peroxidase

4.8 **Wash Buffer** **WASH** **25x** 1 vial
20 ml, concentrated
Dilute content with dist. water to 500 ml total volume.

4.9 **Substrate** **SUB** 1 vial
12 ml TMB solution, ready for use

4.10 **Stop Solution** **STOP** 1 vial
12 ml, ready for use
Contains 0.3 M sulphuric acid

4.11 **Reaction plate** **ACYL-PLATE** 1 piece
for acylation

| | | |
|--|------------------|---------|
| 4.12 Equalizing Reagent | EQUA-REAG | 1 vial |
| lyophilized, dissolve content with 20 ml dist. water, dissolve carefully to minimize foam formation | | |
| 4.13 Adhesive Foil | FOIL | 1 piece |

Additional materials and equipment required but not provided:

- Pipettes (10, 25, 50, 100 and 200 µl)
- Orbital shaker
- Microplate washing device
- Microplate photometer (450 nm)

5. Sample Collection

5.1. Serum and Plasma

The test can be performed with serum as well as with EDTA plasma and heparin plasma. If plasma is to be used care must be taken to get true platelet-free plasma. Otherwise, the Serotonin level has to be related to the number of thrombocytes in the sample. Since the preparation of platelet-free plasma requires special precautions, it is generally recommended to use serum instead of plasma.

Hemolytic and lipemic samples should not be used.

The samples can be stored up to 24 hours at 2 - 8 °C. For a longer storage (up to 6 months) the samples must be frozen at -20 °C

Repeated freezing and thawing should be avoided.

5.2. Urine

The total volume of urine excreted during a 24-hours period should be collected and mixed in a single bottle containing 10 - 15 ml of 6 M hydrochloric acid as preservative. Avoid exposure to direct sun light. Determine the total volume and take an aliquot for the measurement. For patients with suspected kidney disorders the creatinine concentration should be tested, too. Urine samples can be stored at -20 °C for at least 6 months.

6. Preparation of Reagents and Samples

6.1. Microtiter strips **STRIPS**

Before opening the packet of strip wells, allow it to stand at room temperature for at least 10 minutes. After opening, keep any unused wells in the original foil packet with the desiccant provided. Reseal carefully and store at 2-8 °C.

6.2 Wash Buffer **WASH**

Dilute the content with dist. water to a total volume of 500 ml.

For further use the diluted wash buffer must be stored at 2 - 8 °C for a maximum period of 4 weeks.

6.3. Equalizing Reagent **EQUA-REAG**

Dissolve the content with 20 ml dist. water, mix shortly and leave on a roll mixer for 30 minutes. Handle carefully in order to minimize foam formation. The reconstituted Equalizing Reagent should be stored frozen at -20 °C and is stable for a minimum of 1 year.

All other reagents are ready for use.

6.4. Preparation of Samples (Acylation)

The wells of the reaction plate for the acylation can be used only once. So please mark the respective wells before using.

1. Pipette each **10 µl standard 1 - 6**, each **10 µl control 1 & 2**, each **10 µl serum**, **25 µl plasma**, **10 µl urine** or **50 µl cerebrospinal fluid** into the respective wells of the reaction plate.
2. Pipette each **25 µl Acylation Buffer** into all wells.
3. Pipette each **200 µl Equalizing Reagent** into all wells.
4. Pipette each **25 µl Acylation Reagent** into all wells, mix **immediately**.
5. Incubate for 30 minutes at room temperature (approx. 20 °C) on an orbital shaker.
6. Take each 20 µl for the ELISA.

7. Test Procedure ELISA

7.1 Sample Incubation

Pipette each **20 µl prepared Standards 1 to 6**, **20 µl prepared controls** and **20 µl prepared samples** into the respective wells of the coated microtiter strips (duplicates are recommended).

Pipette each **50 µl Antiserum** into all wells.

Cover the plate with adhesive foil and incubate for 60 minutes at room temperature on an orbital shaker.

7.2 Washing

Discard or aspirate the contents of the wells and wash thoroughly with each **250 µl Wash Buffer**. Repeat the washing procedure 3 to 4 times. Remove residual liquid by tapping the inverted plate on clean absorbent paper.

7.3 Conjugate Incubation

Pipette each **100 µl enzyme conjugate** into all wells.

Incubate for 30 minutes at room temperature on an orbital shaker.

7.4 Washing

Repeat step 7.2.

7.5 Substrate Incubation

Pipette each **100 µl Substrate** into all wells and incubate for 20 to 30 minutes at room temperature on an orbital shaker.

7.6 Stopping

Pipette each **100 µl Stop Solution** into all wells.

7.7 Reading

Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer.

8. Calculation of the Results


On a semilogarithmic graph paper the concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Alternatively, the optical density of each standard and sample can be related to the optical density of the zero standard, expressed as the ratio OD/OD_{max}, and then plotted on the y-axis.

The concentration of the controls, urine and serum samples can be read directly from this standard curve by using their average optical density.

The read-off values for plasma samples have to be divided by a factor of 2.4.

The values for cerebrospinal fluid samples have to be divided by a factor of 4.3.

Typical standard curve:



$$y = ((A - D) / (1 + (x / C) ^ B)) + D$$

| | <u>A</u> | <u>B</u> | <u>C</u> | <u>D</u> | <u>R²</u> |
|---|----------|----------|----------|----------|----------------------|
| ○ Std (Standards: Concentration vs MeanValue) | 2,205 | 0,947 | 119,181 | 0,159 | 0,999 |

9. Assay Characteristics

Normal Range

The reference ranges given below should only be taken as a guideline. It is recommended that each laboratory should establish its own normal values.

Plasma (platelet-free): < 10 ng/ml

Serum: Female 80 - 450 ng/ml

Male 40 - 400 ng/ml

Urine: 50 - 250 µg/day

Sensitivity

5 ng/ml

Specificity (Cross Reactivity)

Structural related components were tested for possible interference with the antisera against Serotonin used in the ELISA method. The tested compounds were Tryptamine, Melatonin, 5-HIAA, 5-Methoxytryptamine and L-Tryptophan.

| Substance | ED-50-Value (ng/ml) | Cross Reactivity (%) |
|---------------------|--------------------------------|-----------------------------|
| Serotonin | 131 | 100 |
| Tryptamine | 10,049 | 1.3 |
| 5-Methoxytryptamine | 73,950 | 0.18 |
| Melatonin | > 1,000,000 | < 0.013 |
| 5-HIAA | > 1,000,000 | < 0.013 |
| L-Tryptophan | > 1,000,000 | < 0.013 |

Recovery

Increasing amounts of Serotonin were added to a serum and an urine sample. Each spiked sample was assayed. The analytical recovery of Serotonin was estimated at different concentrations by using the theoretically expected and the actually measured values. The mean recoveries from all concentrations were for serum and urine 107%.

Concentrations in ng/ml

Serum:

| added | measured | expected | % recovery |
|---------------|----------|----------|------------|
| 0 | 47.4 | | |
| 19.6 | 65.3 | 67.0 | 97 |
| 38.5 | 94.4 | 85.9 | 110 |
| 68.4 | 134 | 116 | 116 |
| 96.8 | 160 | 144 | 111 |
| 125 | 205 | 172 | 119 |
| 152 | 220 | 199 | 110 |
| 188 | 235 | 235 | 100 |
| 244 | 322 | 291 | 111 |
| 310 | 418 | 357 | 117 |
| 385 | 456 | 432 | 105 |
| 476 | 481 | 523 | 92 |
| 684 | 729 | 731 | 100 |
| 874 | 932 | 921 | 101 |
| mean recovery | | | 107 |

Urine:

| added | measured | expected | % recovery |
|---------------|----------|----------|------------|
| 0 | 86.4 | | |
| 19.6 | 99.7 | 106 | 94 |
| 38.5 | 137 | 125 | 110 |
| 68.4 | 174 | 155 | 112 |
| 96.8 | 200 | 183 | 109 |
| 125 | 208 | 211 | 98 |
| 152 | 261 | 238 | 109 |
| 188 | 294 | 274 | 107 |
| 244 | 370 | 330 | 112 |
| 310 | 406 | 396 | 102 |
| 385 | 516 | 471 | 109 |
| 476 | 645 | 562 | 115 |
| 684 | 845 | 770 | 110 |
| 874 | 985 | 960 | 103 |
| mean recovery | | | 107 |

Linearity

The linearity of the ELISA method was investigated using different dilutions of a serum and an urine sample. The mean linearity from all dilutions were for serum and urine 103%.

Concentrations in ng/ml

Serum:

| dilution | measured | recalculated value | recovery % |
|----------|----------|--------------------|------------|
| orig. | 1077 | | |
| 4 + 1 | 903 | 862 | 105 |
| 2 + 1 | 721 | 718 | 100 |
| 1 + 1 | 519 | 539 | 96 |
| 1 + 2 | 363 | 359 | 101 |
| 1 + 3 | 296 | 269 | 110 |
| 1 + 4 | 210 | 215 | 97 |
| 1 + 6 | 145 | 154 | 94 |
| 1 + 9 | 117 | 108 | 109 |
| 1 + 14 | 77.7 | 71.8 | 108 |
| 1 + 19 | 57.5 | 53.9 | 107 |

mean linearity

103

Urine:

| dilution | measured | recalculated value | recovery % |
|----------|----------|--------------------|------------|
| orig. | 1121 | | |
| 4 + 1 | 842 | 897 | 94 |
| 2 + 1 | 812 | 747 | 109 |
| 1 + 1 | 585 | 561 | 104 |
| 1 + 2 | 391 | 374 | 105 |
| 1 + 3 | 277 | 280 | 99 |
| 1 + 4 | 196 | 224 | 87 |
| 1 + 6 | 165 | 160 | 103 |
| 1 + 9 | 123 | 112 | 110 |
| 1 + 14 | 89.7 | 74.7 | 120 |

mean linearity

103

Reproducibility

The reproducibility of the ELISA method was investigated by determining the intra-assay-coefficients of variation (cv) by repeated measurements of different serum and urine samples with different Serotonin concentrations.

Intra-Assay Variation

Serum:

| sample | n = | mean value | sd | cv (%) |
|--------|-----|------------|------|--------|
| 1 | 41 | 133 | 8.0 | 6.0 |
| 2 | 41 | 366 | 17.4 | 4.7 |

Urine:

| sample | n = | mean value | sd | cv (%) |
|--------|-----|------------|------|--------|
| 1 | 40 | 128 | 8.4 | 6.6 |
| 2 | 40 | 308 | 16.6 | 5.4 |

10. Literature

- Harenberg, J., Huhle, G., Giese, Ch., Wang, L., Feuring, M., Song, X., Hofmann, U.
Determination of serotonin release from platelets by enzyme immunoassay in the diagnosis of heparin-induced thrombocytopenia. British Journal of Hematology, **109**, 182-186 (2000)
- Balaskas, E., Bamihas, G., Karamouzis, M., Voyiatzis, G., Tourkantonis, A.
Histamine and Serotonin in uremic pruritus: effect of ondansetron in CAPD-pruritic patients
Nephron, **78:395-402** (1998)
- Stratz, T., Samborski, W., Hrycaj, P., Pap, T., Mackiewicz, S., Mennet, P., Müller, W.
Die Serotoninkonzentration im Serum bei Patienten mit generalisierter Tendomyopathie (Fibromyalgie) und chronischer Polyarthritits
Medizinische Klinik, **88**, 458-462 (1993)
- Manz, B. et al.
Radioimmunoassay of human serum serotonin
J. Clin. Chem. Clin. Biochem. **23**, (1985), 657
- Manz, B., Kosfeld H., Belovsky, O., Grill H.-J., Pollow, K.
Direkte radioimmunologische Bestimmung des Serotonins in normalen und pathologischen Seren
Ärztl. Lab. **32** (1986), 135 - 140
- Attanasio, A. et al.
Diurnal rhythms of N-acetylserotonin and serotonin
J. Pineal res. **3** (1986), 251 - 256
- Manz, B., Raynaud, F. and Pevet, P.
Direct radioimmunoassay of N-acetylserotonin and serotonin in the pineal gland of rat and hamster
In: EPSG Newsletter Suppl. 7, Ares Serono Symposia, I. Colloquium of the european pineal study group (1987) p.92
- Kema, I.P. et al.
Influence of a Serotonin- and Dopamine-Rich Diet on Platelet Serotonin Content and Urinary Excretion of Biogenic Amines and Their Metabolites
Clin. Chem. **38/9** (1992), p.1730 – 1736
- Manz, B. et al.
Radioimmunoassay of human serum serotonin
J. Clin. Chem. Clin. Biochem. **23**, (1985), 657
- Kema, I.P. et al.
Improved Diagnosis of Carcinoid Tumors by Measurement of Platelet Serotonin
Clin. Chem. **38/4** (1992), p. 534 - 540

Pipetting Scheme

Sample Preparation

| | Standards | Control | Serum, Urine | Plasma | Cerebr. Fluid |
|--------------------------------|-----------|---------|--------------|--------|---------------|
| Standard 1 - 6 μl | 10 | | | | |
| Control 1 & 2 μl | | 10 | | | |
| Serum, Urine μl | | | 10 | | |
| Plasma μl | | | | 25 | |
| Cerebr. Fluid μl | | | | | 50 |
| Acyl. Buffer μl | 25 | 25 | 25 | 25 | 25 |
| Equalizing Reag. μl | 200 | 200 | 200 | 200 | 200 |
| Acyl. Reagent μl | 25 | 25 | 25 | 25 | 25 |

30 minutes incubation at room temperature

Pipetting Scheme ELISA

| | Standard | Control | Sample |
|------------------------------|----------|---------|--------|
| Standard 1 - 6 μl | 20 | | |
| Control 1 & 2 μl | | 20 | |
| Sample μl | | | 20 |
| Antiserum μl | 50 | 50 | 50 |

60 minutes incubation at room temperature

3 - 4 x washing

| | | | |
|--------------------------------|-----|-----|-----|
| Enzyme Conjugate μl | 100 | 100 | 100 |
|--------------------------------|-----|-----|-----|

30 minutes incubation at room temperature

3 - 4 x washing

| | | | |
|-------------------------|-----|-----|-----|
| Substrate μl | 100 | 100 | 100 |
|-------------------------|-----|-----|-----|

20 - 30 minutes incubation at room temperature

| | | | |
|-----------------------------|-----|-----|-----|
| Stop Solution μl | 100 | 100 | 100 |
|-----------------------------|-----|-----|-----|

Reading of absorbance at 450 nm