



## Instructions for Use

# ADMA / Arginine ELISA


Enzyme Immunoassay  
for the Quantitative Determination of  
**Endogenous Asymmetric Dimethylarginine (ADMA)**  
**and L-Arginine (ARG) in Serum and EDTA-Plasma**

**RUO**

For Research Use Only  
Not for Use in Diagnostic Procedures

**REF** EA215/96

 96

 2 – 8 °C



DLD Gesellschaft für Diagnostika und medizinische Geräte mbH

Adlerhorst 15 • 22459 Hamburg • Germany

Tel +49 40 555 87 10 • Fax +49 40 555 87 111










Internet: <http://www.dld-diagnostika.de> • E-Mail: [contact@dld-diagnostika.de](mailto:contact@dld-diagnostika.de)



## Inhaltsverzeichnis

1	Introduction and Principle of the Test .....	4
2	Precautions .....	5
3	Changes to Declare .....	5
4	Storage and Stability .....	6
5	Contents of the Kit .....	6
6	Sample Collection.....	8
7	Preparation of Reagents.....	9
8	Test Procedure ELISA.....	10
9	Calculation of the Results.....	13
10	Assay Characteristics.....	15
11	Literature .....	17
	Pipetting Scheme – Sample Preparation .....	19
	Pipetting Scheme - ELISA.....	20

## Symbols

	Research use only		
	Content		Expiry Date
	Lot Number		Store at
	Manufactured by		Sufficient for ... determinations
	Catalogue Number		Consult Instructions for Use

The symbols of the components of the kit are described in chapter 4 Contents of the Kit.

## 1 Introduction and Principle of the Test

The ADMA/Arginine ELISA Kit contains reagents for the quantitative determination of derivatized ADMA (asymmetric dimethyl arginine) and arginine in serum or EDTA plasma. Derivatization takes place during sample preparation. ADMA is quantitatively converted into N-acyl-ADMA by the acylation reagent, and arginine is quantitatively converted into N-acyl-arginine.

The ADMA/Arginine ELISA is a competitive enzyme immunoassay. Antigen bound to the solid phase and free antigen in solution compete for a defined number of antibody binding sites. When the system is in equilibrium, the unbound antigen-antibody complex is removed in a washing step and the correspondingly bound complex is detected using an anti-rabbit IgG peroxidase conjugate and determined via the conversion of tetramethylbenzidine (TMB). The TMB/POD reaction is stopped and measured at 450 nm. The concentration of the antigen-antibody complex bound to the solid phase is inversely proportional to the concentration of the antigen in the sample.

## 2 Precautions

- For research use only. Not for use in diagnostic procedures. For professional use only.
- Before performing the test, the valid instructions for use, as included in this kit, should be read completely and the content understood.
- Material of animal origin used in the preparation of the kit have been obtained from certified healthy animals; however, the materials should be handled as potentially infectious.
- Individual components of different lots and kits should not be interchanged. The expiry dates and storage conditions stated on the packaging and the labels of the individual components must be observed.
- When handling the reagents, controls and samples, the current laboratory safety guidelines and good laboratory practice should be observed.
- Wear protective clothing, disposable gloves, and safety goggles while performing the test.
- Some of the components of this test kit contain hazardous substances. These components bear the appropriate hazard symbol on their label. Further information can be found in 4. *Contents of the kit* and on the relevant safety data sheets.
- Avoid any actions that could result in ingestion, inhalation or injection of the reagents. Never pipette by mouth.
- Avoid contact with individual reagents, as these can cause irritation and chemical burns.
- Dispose of waste according to state and local environmental protection regulations.
- Broken glass can cause injury. Be cautious with glass vials.

## 3 Changes to Declare

Version \_8: IFU was reformatted. Additions and changes are highlighted in gray.



Version \_7: Additions and changes are highlighted in gray.





#### 4 Storage and Stability

The kit is shipped at ambient temperature. Upon arrival, store the kit at 2 - 8 °C to keep it stable until its expiry date. Once opened, the kit is stable until its expiry date. The shelf life of the ready-to-use reagents is indicated on the respective bottle label. For stability of prepared reagents refer to 7.

Bring all reagents to room temperature before use and refrigerate immediately after use.

#### 5 Contents of the Kit

<b>MT-Strips</b> 8 wells each, break apart precoated with ADMA; blue colored	<b>STRIPS-ADMA</b>	12 strips
<b>MT-Strips</b> 8 wells each, break apart precoated with L-Arginine; yellow colored	<b>STRIPS-ARG</b>	12 strips
<b>Standards 1 - 6</b> 4 ml each, ready for use	<b>CAL 1 – 6</b>	6 vials
<b>Control 1 &amp; 2</b> 4 ml each, ready for use Range: see QC certificate	<b>CON 1 &amp; 2</b>	2 vials
<b>Acylation Buffer</b> 3.5 ml, ready for use, color coded blue	<b>ACYL-BUFF</b>	1 vial  Warning
<b>Acylation Reagent</b> lyophilized, (see 7.3)	<b>ACYL-REAG</b>	3 vials
<b>Antiserum ADMA</b> 7 ml, ready for use, color coded blue Rabbit-anti-N-acyl-1-ADMA	<b>AS-ADMA</b>	1 vial  Warning

<b>Antiserum Arginine</b> 7 ml, ready for use, color coded yellow Rabbit-anti-N-acyl-1-Arginine	<b>AS-ARG</b>	1 vial	 Warning
<b>Enzyme Conjugate</b> 13 ml, ready for use Anti-rabbit-IgG-peroxidase	<b>CONJ</b>	2 vials	 Warning
<b>Wash Buffer</b> 20 ml, concentrated (50x) (see 7.1)	<b>WASH</b>	2 bottles	
<b>Substrate</b> 13 ml TMB solution, ready for use	<b>SUB</b>	2 vials	
<b>Stop Solution</b> 13 ml, ready for use contains 0.3 M sulphuric acid, not corrosive	<b>STOP</b>	2 vials	
<b>Reaction Plate</b> For acylation	<b>ACYL-PLATE</b>	1 piece	
<b>Equalizing Reagent</b> lyophilized, (see 7.2)	<b>EQUA-REAG</b>	1 vial	
<b>Solvent</b> 10 ml, ready for use, contains DMSO (please note that Solvent reacts with many plastic materials including plastic trays; Solvent does not react with normal pipette tips and with glass devices)	<b>SOLVENT</b>	1 vial	 Warning  Danger
<b>Self-adhesive Foil</b> ready for use	<b>FOIL</b>	4 pieces	

Additional materials and equipment required but not provided:

- Pipettes (10, 20, 25, 50, 100 and 200  $\mu$ l)
- Multichannel pipette or microplate washing device
- Orbital shaker
- Multipipette
- Microplate washing device
- Microplate photometer (450 nm)
- Vortex mixer
- Roll mixer
- Paper towels, pipette tips, timer

## **6 Sample Collection**

The test can be performed with EDTA plasma and serum.

Haemolytic and lipemic samples should not be used.

The samples can be stored up to 6 hours at 2 – 8 °C. For a longer storage (up to 18 months) the samples must be kept frozen at -20 °C

Repeated freezing and thawing should be avoided.

## 7 Preparation of Reagents

Equilibrate reagents to room temperature.

### 7.1 Wash Buffer

Dilute the content (20 ml) of 50x concentrated Wash Buffer **WASH** with dist. water to a total volume of 1,000 ml, mix briefly. For further use, the diluted wash buffer must be stored at 2 – 8 °C for a maximum period of 4 weeks.

Should the kit be used in several runs, then prepare only the required amount of wash buffer for each run.

### 7.2 Equalizing Reagent

Dissolve the content of **EQUA-REAG** with 21 ml dist. water, mix shortly and leave on a roll mixer or orbital shaker for 20 minutes. Handle carefully in order to minimize foam formation. The reconstituted Equalizing Reagent should be stored frozen at -20 °C and is stable until expiry date of the kit.

### 7.3 Acylation Reagent

Remove the required amount of vials of Acylation Reagent **ACYL-REAG** from the foil pouch, leave the remaining vials inside together with the desiccant and close the pouch carefully. Reconstitute each vial of lyophilized Acylation Reagent with 3 mL of Solvent **SOLVENT** and shake for minimum 10 minutes on an orbital shaker. The Acylation Reagent should be freshly prepared immediately before performing the test and is then stable for approx. 3 hours. The second and third vial allow a second and third run, respectively, of the test. If the whole kit is to be used in one run, it is recommended to pool the dissolved contents of two vials of Acylation Reagent. Discard the remaining reconstituted reagent after use.

**All other reagents are ready for use.**

## 8 Test Procedure ELISA

Equilibrate reagents to room temperature.

Duplicates are recommended.

The wells of the reaction plate **ACYL-PLATE** for the acylation can be used only once. Please mark the respective wells before use to avoid repeated use.

### 8.1 Preparation of Samples (Acylation)

1. Pipette each 20 µl Standard 1 to 6 **CAL 1-6**, Control 1 & 2 **CON 1 & 2** and sample into the respective wells of the Reaction Plate **ACYL-PLATE**.
2. Pipette 20 µl Acylation Buffer **ACYL-BUFF** into each well.
3. Pipette 200 µl Equalizing Reagent **EQUA-REAG** (see 7.2) into each well. Mix for 10 seconds on an orbital shaker at medium speed.
4. Prepare Acylation Reagent **ACYL-REAG** (see 7.3) just before use. For pipetting please use a multipipette or similar device, fill the syringe directly from the vial with dissolved Acylation Reagent and add well by well. Please note that dissolved Acylation Reagent reacts with many plastic materials including plastic trays. It does not react with normal pipette tips and with glass devices. Pipette 50 µl prepared Acylation Reagent **ACYL-REAG** into each well and continue with the next step, immediately. The color changes to violet.
5. Incubate for 20 minutes at room temperature (20 - 25 °C) on an orbital shaker at medium speed.

Take 25 µl each for the ADMA ELISA and 10 µl each for the Arginine ELISA.

## 8.2 ADMA ELISA

1. Pipette each 25 µl prepared Standards, prepared Controls and prepared samples into the respective wells of the coated microtiter strips **STRIPS-ADMA** (blue colored).  
Leave remaining microtiter strips in the foil pouch together with the desiccant and close thoroughly.
2. Pipette 50 µl ADMA-Antiserum **AS-ADMA** into each well.
3. Cover the plate with adhesive foil **FOIL** and incubate for 90 minutes at room temperature (20 - 25 °C) on an orbital shaker at medium speed.
4. Discard or aspirate the contents of the wells and wash thoroughly with 300 µl Wash Buffer **WASH** (see 7.1) per well. Discard or aspirate the contents of the wells and remove residual liquid by tapping the inverted plate on a clean absorbent paper towel. Repeat the washing procedure 3 times.  
Alternatively, a washing device may be used.
5. Pipette 100 µl Enzyme Conjugate **CONJ** into each well.
6. Incubate for 30 minutes at room temperature (20 – 25 °C) on an orbital shaker at medium speed.
7. Washing: Repeat step 4.
8. Pipette 100 µl Substrate **SUB** into each well.
9. Incubate for 25 ± 5 minutes at room temperature (20 – 25 °C) on an orbital shaker at medium speed.
10. Pipette 100 µl Stop Solution **STOP** into each well and shake on an orbital shaker for minimum 10 seconds at medium speed.
11. Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer within 15 minutes.

### 8.3 Arginine-ELISA

1. Pipette each 10 µl prepared Standards, prepared controls and prepared samples into the respective wells of the coated microtiter strips **STRIPS-ARG** (yellow colored).  
Leave remaining microtiter strips in the foil pouch together with the desiccant and close thoroughly.
2. Pipette each 50 µl Arginine-Antiserum **AS-ARG** into each well.
3. Cover the plate with adhesive foil **FOIL** and incubate for 90 minutes at room temperature (20 - 25 °C) on an orbital shaker at medium speed.
4. Discard or aspirate the contents of the wells and wash thoroughly with 300 µl Wash Buffer **WASH** (see 7.1) per well. Discard or aspirate the contents of the wells and remove residual liquid by tapping the inverted plate on a clean absorbent paper towel. Repeat the washing procedure 3 times.  
Alternatively, a washing device may be used.
5. Pipette 100 µl Enzyme Conjugate **CONJ** into each well.
6. Incubate for 30 minutes at room temperature (20 – 25 °C) on an orbital shaker at medium speed.
7. Washing: Repeat step 4.
8. Pipette 100 µl Substrate **SUB** into each well.
9. Incubate for 25 ± 5 minutes at room temperature (20 – 25 °C) on an orbital shaker at medium speed.
10. Pipette 100 µl Stop Solution **STOP** into each well and shake on an orbital shaker for minimum 10 seconds at medium speed.
11. Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer within 15 minutes.

## 9 Calculation of the Results

Standard	1	2	3	4	5	6
ADMA (µmol/l)	0	0.2	0.45	0.7	1	3
Arginine (µmol/l)	5	15	35	70	120	300

The concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Cubic spline, 4 parameter or similar iteration procedures are recommended for evaluation of the standard curve.

The concentration of the controls and samples in µmol/l can be read directly from this standard curve by using their average optical density.

Serum arginine levels may be approximately 50% higher than the corresponding EDTA plasma levels in the same donor.


Conversion:

ADMA: 1 µmol/l = 202 ng/ml


Arginine: 1 µmol/l = 174 ng/ml

Typical standard curve (do not use for calculation of results):

### ADMA ELISA



## Arginine ELISA:



$y = ( (A - D) / (1 + (x/C)^B) ) + D$ :    A    B    C    D    R<sup>2</sup>  
 ○ Std (Standards: Concentration vs MeanValue)    2,248    0,975    32,2    -0,071    1

Quality Control: Test results are valid only if the kit controls are within the ranges specified on the QC Certificate. Otherwise, the test should be repeated.

## 10 Assay Characteristics

### 10.1 Sensitivity

	Lower Limit of Detection	Calculation
ADMA	0.03 $\mu\text{mol} / \text{l}$	$\text{OD}_{\text{Cal1}} - 3 \times \text{SD}$
Arginin	6 $\mu\text{mol} / \text{l}$	$\text{OD}_{\text{Cal1}} - 3 \times \text{SD}$

### 10.2 Specificity (Cross Reactions)

#### ADMA

Substance	Cross Reactivity (%)
ADMA	100
SDMA	0.05
Monomethylarginine (NMMA)	1.93
Homoarginine	< 0.01
Arginine	0.03

#### Arginine

Substance	Cross Reactivity (%)
Arginine	100
ADMA	< 0.37
Homoarginine	2.92
SDMA	0.88

### 10.3 Recovery

#### ADMA

	Range ( $\mu\text{mol} / \text{l}$ )	Mean (%)	Range (%)
EDTA-Plasma	0.43 – 1.55	99	90 - 107
Serum	0.54 – 1.72	92	87 - 102

#### Arginine

	Range ( $\mu\text{mol} / \text{l}$ )	Mean (%)	Range (%)
EDTA-Plasma	48 – 163	97	93 - 100
Serum	82 – 211	100	96 - 103

### 10.4 Linearity

#### ADMA

	Range ( $\mu\text{mol} / \text{l}$ )	Highest Dil.	Mean (%)	Range (%)
EDTA-Plasma	0.23 – 1.53	1 : 6 with water	99	92 - 105

**Arginine**

	Range ( $\mu\text{mol} / \text{l}$ )	Highest Dil.	Mean (%)	Range (%)
EDTA-Plasma	28 – 193	1 : 6 with water	102	94 - 106

**10.5 Reproducibility****ADMA**

	Range ( $\mu\text{mol} / \text{l}$ )	Intra-Assay-CV
EDTA-Plasma	0.58 – 1.04	4.9 – 5.4 %

	Range ( $\mu\text{mol} / \text{l}$ )	Inter-Assay-CV
EDTA-Plasma	0.57 – 1.34	4.3 – 9.6 %

**Arginine**

	Range ( $\mu\text{mol} / \text{l}$ )	Intra-Assay-Vk
EDTA-Plasma	56 – 125	3.6 – 2.3 %

	Range ( $\mu\text{mol} / \text{l}$ )	Inter-Assay-Vk
EDTA-Plasma	53 – 170	3.2 – 6.3 %

**10.6 Method Comparison****ADMA**

	Method	Correlation
Serum + EDTA-Plasma	LC/MS	$Y = 0.99 \times \text{LC/MS} + 0.02$ ; $R = 0.983$ ; $N = 32$

**Arginine**

	Method	Correlation
Serum + EDTA-Plasma	LC/MS	$Y = 0.95 \times \text{LC/MS} - 0.68$ ; $R = 0.991$ ; $N = 32$

**10.7 Calibration**

The calibration is carried out by weighing the pure substance.

**10.8 Limitations of Method**

The results are to be used for research use only.

Samples measured above the highest standard must be diluted with distilled water and reassayed (see 10.4) The values of diluted samples must be multiplied by the appropriate dilution factor.

**10.9 Interferences**

Hemolytic, lipemic and icteric specimens should not be used.

## 11 Literature

### 11.1 Literature using the ADMA-ELISA from DLD Diagnostika

Schulze F, Wesemann R, Schwedhelm E, Sydow K, Albsmeier J, Cooke JP, Böger RH.

**Determination of ADMA using a novel ELISA assay.**

Clin. Chem. Lab. Med. 2004; 42: 1377-1383

Krempl TK, Kähler J, Maas R, Silberhorn L, Meinertz T, Böger RH.

**Elevation of asymmetric dimethylarginine (ADMA) in patients with unstable angina and recurrent cardiovascular events.**

Eur. Heart J. 2005; 26: 1846-1851

Schulze F, Maas R, Freese R, Schwedhelm E, Silberhorn L, Böger RH.

**Determination of a reference value for N,N-dimethyl-L-arginine in 500 subjects.**

Eur. J. Clin. Invest. 2005; 35 : 622-626

Schnabel R, Blankenberg S, Lubos E, Lackner KJ, Rupprecht HJ, Espinola-Klein C, Jachmann N, Post F, Peetz D, Bickel C, Cambien F, Tiret L, Münzel T.

**Asymmetric dimethylarginine and the risk of cardiovascular events and death in patients with coronary artery disease: results from the AtheroGene Study.**

Circ. Res. 2005; 97: e53-59

O'Dwyer MJ, Dempsey F, Crowley V, Kelleher D, McManus R, Ryan T.

**Septic shock correlates with ADMA levels which may be influenced by a polymorphism in DDAH II: a prospective observational study.**

Crit. Care 2006; 10: (5): R139

Antoniades C, Tousoulis D, Marinou K, Vasiliadou C, Tentolouris C, Bouras G, Pitsavos C, Stefanidis C.

**Asymmetrical dimethylarginine regulates endothelial function in methionine-induced but not in chronic homocystinemia in humans: effect of oxidative stress and proinflammatory cytokines**

Am. J. Clin. Nutr. 2006; 84: 781-788

Wang TZ., Chen WJ., Cheng WC., Lin JW., Chen MF., Lee YT.

**Relation of improvement in endothelium-dependent flowmediated vasodilation after Rosiglitazone to changes in asymmetric dimethylarginine, endothelin-1, and C-reactive protein in nondiabetic patients with the metabolic syndrome**

Am. J. Cardiol. 2006; 9: 1057-1062

Wanby P., Nilsson I., Brudin L., Nyhammar I., Gustafsson I., Carlsson M.

**Increased plasma levels of asymmetric dimethylarginine in patients with carotid stenosis: no evidence for the role of the common FABBP2 A54T gene polymorphism**

Acta Neurol. Scand. 2007; 115: 90-96

Konishi H, Sydow K, Cooke JP.

**Dimethylarginine dimethylaminohydrolase promotes endothelial repair after vascular injury**

J. Am. Coll. Cardiol. 2007; 49: 1099-1105

Iribarren C, Husson G, Sydow K, Wang BY, Sidney S, Cooke JP.

**Asymmetric dimethyl-arginine and coronary artery calcification in young adults entering middle age: the CARDIA Study**

Eur. J. Cardiovasc. Prev. Rehabil. 2007; 14:222-229

Melikian N, Wheatcroft SB, Ogah OS, Murphy C, Chowienczyk PJ, Wierzbicki AS, Sanders TA, Jiang B, Duncan ER, Shah AM, Kearney MT.

**Asymmetric dimethylarginine and reduced nitric oxide bioavailability in young Black African men**  
Hypertension 2007; 49: 873-877

Horowitz JD, Heresztyn T.

**An overview of plasma concentrations of asymmetric dimethylarginine (ADMA) in health and disease and in clinical studies: Methodological considerations.**

J. Chromatogr. B Analyt. Technol. Biomed. Life Sci. 2007; epub ahead of print

Korish AA, Arafah MM.

**Catechin combined with vitamins C and E ameliorates insulin resistance (IR) and atherosclerotic changes in aged rats with chronic renal failure (CRF)**

Arch. Gerontol. Geriatr. 2007; in press

Charitidou C, Farmakiotis D, Zournatzi V, Pidonia I, Pegiou T, Karamanis N, Hatzistilianou M, Katsikis I, Panidis D.

**The administration of estrogens, combined with anti-androgens, has beneficial effects on the hormonal features and asymmetric dimethyl-arginine levels, in women with the polycystic ovary syndrome**

Atherosclerosis 2007; in press

## 11.2 General Literature

Vallance P, Leone A, Calver A, Collier J, Moncada S.

**Accumulation of an endogenous inhibitor of NO synthesis in chronic renal failure**

Lancet 1992; 339: 572 - 575

Stühlinger M, Abbasi F, Chu JW, Lamendola C, McLaughlin TL, Cooke JP, Reaven GM, Tsao PS.

**Relationship between insulin resistance and an endogenous nitric oxide synthase inhibitor**

J. Am. Med. Assoc. 2002; 287: 1420-1426

Zoccali C, Bode-Böger SM, Mallamaci F, Benedetto FA, Tripepi G, Malatino L, Cataliotti A, Bellanuova I, Fermo I, Frölich JC, Böger RH.

**Asymmetric dimethylarginine (ADMA): An endogenous inhibitor of nitric oxide synthase predicts mortality in end-stage renal disease (ESRD)**

Lancet 2001; 358: 2113-2117

Nijveldt RJ, Teerlink T, Van der Hoven B, Siroen MP, Kuik DJ, Rauwerda JA, van Leeuwen PA.

**Asymmetrical dimethylarginine (ADMA) in critically ill patients: high plasma ADMA concentration is an independent risk factor of ICU mortality**

Clin. Nutr. 2003; 22: 23-30

Savvidou MD, Hingorani AD, Tsikas D, Frolich JC, Vallance P, Nicolaidis KH.

**Endothelial dysfunction and raised plasma concentrations of asymmetric dimethylarginine in pregnant women who subsequently develop pre-eclampsia**

Lancet 2003; 361: 1511-1517

Böger RH.

**The emerging role of asymmetric dimethylarginine as a novel cardiovascular risk factor**

Cardiovasc. Res. 2003; 59: 824-833

Lu TM, Ding YA, Lin SJ, Lee WS, Tai HC.

**Plasma levels of asymmetrical dimethylarginine and adverse cardiovascular events after percutaneous coronary intervention.**

Eur Heart J. 2003; 24: 1912-1919

### Pipetting Scheme – Sample Preparation

For ADMA and Arginine together in ACYL-PLATE

		Standards	Controls	Plasma	Serum
ACYL-PLATE:					
CAL 1 - 6	μl	20			
CON 1 & 2	μl		20		
EDTA-Plasma	μl			20	
Serum					20
ACYL-BUFF	μl	20	20	20	20
EQUA-REAG	μl	200	200	200	200

Shake plate for 10 seconds

ACYL-REAG	μl	50	50	50	50
-----------	----	----	----	----	----

**Immediately**, shake for 20 minutes at room temperature

Take **25 μl** each for ADMA ELISA

Take **10 μl** each Arginin ELISA

### Pipetting Scheme - ELISA

For ADMA and Arginine in separate microtiter plates:

		ADMA (blue) STRIPS-ADMA			Arginin (yellow) STRIPS-ARG		
		Acyl. Stand.	Acyl. Ctrl.	Acyl. Samples	Acyl. Stand.	Acyl. Ctrl.	Acyl. Samples
Transfer from ACYL-PLATE into STRIPS		25	25	25	10	10	10
AS-ADMA (blue)	µl	50	50	50	-	-	-
AS-ARG (yellow)	µl	-	-	-	50	50	50

Seal STRIPS with FOIL

Shake for 90 minutes at room temperature

Wash 4 x with 300 µl WASH per well

CONJ	µl	100	100	100	100	100	100
------	----	-----	-----	-----	-----	-----	-----

Shake for 30 minutes at room temperature

Wash 4 x with 300 µl WASH per well

SUB	µl	100	100	100	100	100	100
-----	----	-----	-----	-----	-----	-----	-----

Shake for 25 ± 5 minutes at room temperature

STOP	µl	100	100	100	100	100	100
------	----	-----	-----	-----	-----	-----	-----

Shake plate for minimum 10 seconds

Read absorbance at 450 nm (ref. 570 – 650 nm)