



## Instruction for use

# Histamine ELISA

Enzyme Immunoassay  
for the Quantitative Determination of  
Histamine in Plasma, Urine and Cell Culture Media

IVD

CE

REF

EA213/96



12 x 8



2 – 8 °C



DLD Gesellschaft für Diagnostika und medizinische Geräte mbH  
Adlerhorst 15 • 22459 Hamburg • Germany

Tel +49-40-555 87 10 • Fax +49-40-555 87 111











Internet: <http://www.dld-diagnostika.de> • E-Mail: [contact@dld-diagnostika.de](mailto:contact@dld-diagnostika.de)



## Contents

1. Introduction and Principle of the Test	Page 5
2. Precautions	Page 6
3. Storage and Stability	Page 6
4. Contents of the Kit	Page 6
5. Specimen Collection and Storage	Page 8
6. Preparation of Samples (Acylation)	Page 9
7. Test Procedure ELISA	Page 11
8. Calculation of the Results	Page 12
9. Assay Characteristics	Page 13
10. Literature	Page 15
Pipetting Scheme	Page 16

## Symbols

	In Vitro Diagnostic Medical Device		EC Declaration of Conformity
	Contents		Expiry Date
	Lot Number		Store
	Manufactured by		Sufficient for
	Catalogue Number		Consult Instructions

## Hazard Pictograms





## 1. Introduction and Principle of the Test

Histamine ( $\beta$ -imidazole-ethylamine) a biogenic amine, is a product of the histidine metabolism. It is produced by decarboxylation of histidine.

Histamine is widely distributed in mammalian tissues. It's bound to heparin (as inactive form) and stored in the granules of basophilic leukocytes and mast cells and is actively released as required. These cells, if sensitized by IgE antibodies attached to their membranes, degranulate when exposed to the appropriate antigen.

Histamine plays a major role in the initial phase of an anaphylactic reaction.

The quantification of histamine in plasma after allergen administration is of clinical interest.

Histamine is part of the immune response to foreign pathogens and it increases the permeability of the capillaries to white blood cells and other proteins, in order to allow them to engage foreign invaders in the affected tissues. Responsible for the biological effects of histamine in tissue are the activation of different surface receptors, for instance H1, H2 and H3.

Histamine is involved in the regulating physiological function in the gut and acting as a neurotransmitter.

The competitive Histamine ELISA kit uses the microtitre plate format. Histamine is bound to the solid phase of the microtiter plate. Acylated histamine and solid phase bound histamine compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. The antibody bound to the solid phase histamine is detected by anti-rabbit/peroxidase. The substrate TMB / peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase histamine is inversely proportional to the histamine concentration of the sample.

## 2. Precautions

- For in vitro diagnostic use only.
- Disposable gloves and safety glasses should be used.
- All reagents of human origin used in this kit are tested for HIV I/II antibodies, HCV and HBsAg and found to be negative. However, because no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as potentially biohazardous materials.
- Material of animal origin used in the preparation of the kit have been obtained from certified healthy animals but these materials should be handled as potentially infectious.
- Some components of this kit are containing hazardous reagents. These components are marked with the adequate hazard label.

## 3. Storage and Stability

On arrival, store the kit at 2-8 °C. Once opened the kit is stable until its expiry date. For stability of prepared reagents refer to Preparation of Reagents. Do not use components beyond the expiration date shown on the kit labels. Do not mix various lots of any kit component within an individual assay.




## 4. Contents of the Kit

- 4.1 **MT-Strips** **STRIPS** 12 strips  
8 wells each, break apart  
precoated with histamine

- 4.2 **Standards 1 - 6** **CAL 1-6** 6 vials  
Each 4 ml, ready for use  
Concentrations:

Standard	1	2	3	4	5	6
ng/ml	0	0.2	0.6	2	6	25

- 4.3 **Control 1 & 2** **CON 1 & 2** 2 vials  
Each 4 ml, ready for use  
Range: see q.c. certificate

4.4	<b>Acylation Buffer</b> 6 ml, colour coded blue, ready for use	<b>ACYL-BUFF</b>	1 vial
4.5	<b>Acylation Reagent</b> lyophilised, dissolve content in 1.5 ml <b>Solvent</b>	<b>ACYL-REAG</b>	3 vials
4.6	<b>Solvent</b> 5.5 ml solvent to dissolve the Acylation reagent Contains acetone, ready for use	<b>SOLVENT</b>	1 vial
	 Warning  Danger		
4.7	<b>Antiserum</b> 5.5 ml, ready for use, colour coded yellow rabbit-anti-N-acyl-histamine	<b>AS</b>	1 vial
4.8	<b>Enzyme Conjugate</b> 12 ml, ready for use goat anti-rabbit-IgG-peroxidase	<b>CONJ</b>	1 vial
	 Warning		
4.9	<b>Wash Buffer</b> 20 ml, 50x concentrated Dilute content with distilled water to 1 litre total volume	<b>WASH</b>	1 vial
4.10	<b>Substrate</b> 12 ml TMB solution, ready for use	<b>SUB</b>	1 vial
4.11	<b>Stop Solution</b> 12 ml, ready for use contains 0.3 M sulphuric acid	<b>STOP</b>	1 vial
4.12	<b>Reaction Plate</b> for acylation	<b>ACYL-PLATE</b>	2 plates
4.13	<b>Equalizing Reagent</b> lyophilized, dissolve contents with distilled water, volume: see vial label	<b>EQUA-REAG</b>	1 vial

Additional materials and equipment required but not provided:

- Pipettes 20, 30, 50 and 100 µl
- Orbital shaker
- Multichannel pipette or Microplate washing device
- Microplate photometer (450 nm)
- Distilled water

## **5. Specimen Collection and Storage**

The test can be performed with EDTA or Heparin plasma, urine and cell culture media.

Repeated freezing and thawing of samples should be avoided.

### **Plasma**

EDTA or Heparin plasma can be used. Haemolytic and lipaemic samples should not be used.

The samples can be stored up to 6 hours at 2 - 8 °C. For a longer storage (up to 6 months) the samples must be frozen at -20 °C

### **Urine**

Spontaneous urine can be used for this test as well as collected urine. In this case the total volume of urine excreted during a 24-hours period should be collected and mixed in a single bottle containing 10 - 15 ml of 6 M hydrochloric acid as preservative. Avoid exposure to direct sun light. Determine the total volume and take an aliquot for the measurement. For patients with suspected kidney disorders the kreatinine concentration should be tested, too. Urine samples can be stored at -20 °C for at least 6 months.

Urine samples have to be diluted 1:15 with dist. water before assay.

### **Cell Culture Media**

Media like DMEM and RPMI can be used in the test. Other media have to be checked by the user.



## 6. Preparation of Reagents and Samples

### 6.1. Preparation of Reagents

#### Wash Buffer

#### WASH

Dilute the content with dist. water to a total volume of 1,000 ml.

The diluted wash buffer has to be stored at 2 - 8 °C for a maximum period of 4 weeks. For longer storage freeze at -20 °C.

#### Equalizing Reagent

#### EQUA-REAG

Dissolve the contents with dist. water (for volume refer to vial label), mix shortly and leave on a roll mixer for minimum 20 minutes. Handle with care in order to minimize foam formation. The reconstituted Equalizing Reagent should be stored frozen at -20 °C and is stable until expiry date printed on vial label.

#### Acylation Reagent

#### ACYL-REAG

Dissolve the content of one bottle in 1.5 ml **Solvent** and shake for 5 minutes on an orbital shaker. The Acylation Reagent has always to be prepared immediately before use. After use the reagent has to be discarded.

The second and third vial allows a second and third run of the test, respectively. If the whole kit is to be used in one run it is recommended to pool the dissolved contents of the three vials of Acylation Reagent.

Please note that solvent reacts with many plastic materials including plastic trays; solvent does not react with normal pipette tips and with glass devices

Solvent is volatile and the dissolved Acylation Reagent evaporates quickly. Therefore, please do not use a tray with big surface together with a multichannel pipette for pipetting Acylation Reagent. Rather, use an Eppendorf multipipette (or similar device), fill the syringe directly from the vial with dissolved Acylation Reagent and add well by well.

All other reagents are ready for use.

## 6.2. Preparation of Samples (Acylation)

Allow reagents and samples to reach room temperature.

Determinations in duplicates are recommended.

The wells of the reaction plate for the acylation can be used only once. So please mark the respective wells before using.

1. Pipette each 50 µl standard 1 - 6, 50 µl control 1 & 2, 50 µl EDTA plasma samples, 20 µl Heparin plasma, 50 µl urine samples (diluted 1:15 with dist. water) and 50 µl cell culture media samples into the respective wells of the reaction plate.
2. Pipette each 50 µl Acylation Buffer into all wells.
3. Pipette 50 µl dist. water in all wells containing plasma samples.
4. Pipette each 50 µl Equalizing Reagent into wells containing standards, controls, urine samples and cell culture media samples. Pipette each 30 µl Equalizing Reagent into wells containing Heparin plasma. Do not pipette into wells containing EDTA plasma samples. Mix the reaction plate for 10 seconds.
5. Pipette each 10 µl Acylation Reagent into all wells and continue with step 6. immediately. Colour changes to violet.

Please note that solvent reacts with many plastic materials including plastic trays; solvent does not react with normal pipette tips and with glass devices

Solvent is volatile and the dissolved Acylation Reagent evaporates quickly. Therefore, please do not use a tray with big surface together with a multichannel pipette for pipetting Acylation Reagent. Rather, use an Eppendorf multipipette (or similar device), fill the syringe directly from the vial with dissolved Acylation Reagent and well by well.

6. Incubate for 15 minutes at room temperature on an orbital shaker with medium frequency. Do not cover the wells or the plate; leave the plate open on the shaker.
7. Pipette 50 µl Antiserum into all wells.
8. Incubate for 30 minutes at room temperature on an orbital shaker with medium frequency. Do not cover the wells or the plate; leave the plate open on the shaker.

Take each 50 µl for the ELISA.

## 7. Test Procedure ELISA

Allow reagents and samples to reach room temperature.  
Determinations in duplicates are recommended.

1. Pipette each 50  $\mu$ l prepared Standards 1 to 6, Controls and Samples into the respective wells of the coated microtiter strips.
2. Incubate for 60 minutes at room temperature (20 – 25 °C) on an orbital shaker with medium frequency.  
Alternative: Shake the plate briefly manually and incubate 90 minutes without shaking.
3. Discard or aspirate the contents of the wells, add each 300  $\mu$ l Wash Buffer, again discard or aspirate the contents of the wells. Remove residual liquid by tapping the inverted plate on clean absorbent paper.  
Repeat the washing procedure 4 times.
4. Pipette each 100  $\mu$ l enzyme conjugate into all wells.
5. Incubate for 20 minutes at room temperature on an orbital shaker with medium frequency.  
Alternative: Shake the plate briefly manually and incubate 25 minutes without shaking
6. Washing: Repeat step 3.
7. Pipette each 100  $\mu$ l Substrate into all wells.
8. Incubate for 20  $\pm$  5 minutes at room temperature (20 - 25 °C) on an orbital shaker with medium frequency.  
Alternative: Shake the plate briefly manually and incubate 20  $\pm$  5 minutes without shaking
9. Pipette each 100  $\mu$ l Stop Solution into all wells.
10. Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer within 15 minutes.

## 8. Calculation of the Results

On a semilogarithmic graph paper the concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Alternatively, the optical density of each standard and sample can be related to the optical density of the zero standard, expressed as the ratio  $OD/OD_{max}$ , and then plotted on the y-axis. Evaluation by 4 parameter iteration or cubic spline is recommended.


The concentration of the controls and plasma samples and cell culture media can be read directly from this standard curve by using their average optical density.

The read concentration of Heparin plasma samples has to be multiplied by a factor of 2.5.

The read concentration of urine samples has to be multiplied by a factor of 15.

Conversion: 1 ng/ml corresponds to 9,0 nmol/l

### Typical standard curve:



## 9. Assay Characteristics

### 9.1 Normal Range

	Reference Range
EDTA Plasma	< 1 ng/ml
Heparin Plasma	< 4.5 ng/ml
Urin	< 45 µg/day

The reference range given below should only be taken as a guideline. It is recommended that each laboratory should establish its own normal values.

### 9.2 Sensitivity

The lower limit of detection was determined by taking the 2fold standard deviation of the absorbance of the Zero Reference and reading the corresponding value from the standard curve.

	Sensitivity
EDTA Plasma	0.06 ng/ml
Heparin Plasma	0.15 ng/ml
Urin	0.9 ng/ml

### 9.3. Specificity (Cross Reactivity)

Structural related components were tested for possible interference with the antisera against histamine used in the ELISA method.

Substanz	Cross Reactivity (%)
histamine	100
1-methyl histamine	0.054
3-methyl histamine	0.13
1-methyl-4-imidazol-acetic acid	< 0.0001
imidazol-4-acetic acid	< 0.0002
L-histidine	< 0.0001

## 9.4. Recovery

Increasing amounts of histamine were added to each sample. Each spiked sample was assayed. The analytical recovery of histamine was estimated at different concentrations by using the theoretically expected and the actually measured values.

Concentrations in ng/ml

	Range (ng/ml)	Mean (%)	Range (%)
<b>EDTA Plasma</b>	0.6 – 13.4	101	93 - 111
<b>Heparin Plasma</b>	0.8 – 36.0	104	87 - 112
<b>Urine</b>	6.1 – 140.6	98	94 - 103
<b>Cell Culture Media</b>	1.0 – 12.9	104	91 - 121

## 9.5. Linearity

The linearity of the ELISA method was investigated using different dilutions of a sample.

Concentrations in ng/ml

	Range(ng/ml)	Highest Dilution	Mean (%)	Range (%)
<b>EDTA Plasma</b>	0.5 – 10.0	1 : 20 Equalizing Reagent	106	96 - 111
<b>Heparin Plasma</b>	0.9 – 13.7	1 : 15 Equalizing Reagent	102	93 - 109
<b>Urine</b>	7 – 142	1 : 20 dist. water	96	81 - 102
<b>Cell Culture Media</b>	1.1 – 10.3	1 : 10 dist. water	106	99 - 111

## 9.6. Reproducibility

The reproducibility of the ELISA method was investigated by measuring the intra- and inter-assay-coefficients of variation (cv).

Concentrations in ng/ml

	Range (ng/ml)	Intra-Assay-cv (%)	Range (ng/ml)	Inter-Assay-cv (%)
<b>EDTA Plasma</b>	1.2 – 8.7	6.1 – 6.5	1.1 – 3.3	6.2 – 7.3
<b>Heparin Plasma</b>	2.5 – 11.8	6.3 – 5.0	2.1 – 10.8	8.9 – 4.4
<b>Urin</b>	24.1 – 89.6	6.6 – 5.7	15.7 – 43.7	7.2 – 11.3
<b>Zellkultur</b>	1.5 – 5.1	6.3 – 8.6	1.3 – 4.1	10.5 – 6.5

## 10. Literature

- Nettis, E.; Colanardi, A.; Ferrannini, A. (2005):  
**Antihistamines as Important Tools for Regulating Inflammation**  
*Curr. Med. Chem. – Anti-Inflammatory & Anti-Allergy Agents*, 4, 81-89
- Matsumoto, J.; Matsuda, H. (2002):  
**Mast-cell-dependent histamine release after praziquantel treatment of *Schistosoma japonicum* infection: implications for chemotherapy-related adverse effects**  
*Parasitol Res* 88: 888–893
- Belic, A.; Grabnar, I.; Karba, R.; et al. (1999):  
**Interdependence of histamine and methylhistamine kinetics: modelling and simulation approach**  
*Computers in Biology and Medicine* 29, 361-375
- Martens-Lobenhoffer, J.; Neumann, H. (1999):  
**Determination of 1-methylhistamine and 1-methylimidazoleacetic acid in human urine as a tool for the diagnosis of mastocytosis**  
*Journal of Chromatography B*, 721, 135–140
- Prell, G.; Green, J.; Elkashef, A. (1996):  
**The relationship between urine excretion and biogenic amines and their metabolites in cerebrospinal fluid of schizophrenic patients**  
*Schizophrenia Research* 19, 171-176
- Eberlein-König, B.; Ullmann, S.; Thomas, P.; et al. (1995):  
**Tryptase and histamine release due to a sting challenge in bee venom allergic patients treated successfully or unsuccessfully with hyposensitization**  
*Clinical and Experimental Allergy*, Volume 25, pages 704-712
- Koller, D.; Rosenkranz, A.; Pirker, C.; et al. (1992):  
**Assessment of histamine release from basophils in whole blood by benzylpenicilloyl poly-L-lysine in penicillin-sensitized patients**  
*Allergy*: 47: 459-462.
- Marquardt, D.; Wasserman, S. (1982):  
**Mast Cells in Allergic Diseases and Mastocytosis**  
*West J Med*; 137:195-212
- Butchers, P.; Vardey, C.; Skidmore, I.; et al. (1980):  
**Histamine-Containing Cells from Bronchial Lavage of Macaque Monkeys. Time Course and Inhibition of Anaphylactic Histamine Release**  
*Int. Archs Allergy appl. Immun.* 62: 205-212

## Pipetting Scheme Sample Preparation

		Standards	Control	EDTA Plasma	Heparin Plasma	Urine (dil.)	Medium
Standard 1 - 6	µl	50					
Control 1 & 2	µl		50				
EDTA Plasma	µl			50			
Heparin Plasma	µl				20		
Urine (1:15 dil.)	µl					50	
Medium	µl						50
Acyl. Buffer	µl	50	50	50	50	50	50
Dist. Water	µl			50	50		
Equalizing Reag.	µl	50	50		30	50	50

shake for 10 seconds

Acyl. Reagent	µl	10	10	10	10	10	10
---------------	----	----	----	----	----	----	----

**Immediately** shake 15 minutes at room temperature, leave plate open

Antiserum	µl	50	50	50	50	50	50
-----------	----	----	----	----	----	----	----

shake 30 minutes at room temperature, leave plate open

take each 50 µl for the ELISA

## Pipetting Scheme ELISA

		Standards	Control	Sample
Standard 1 - 6	µl	50		
Control 1 & 2	µl		50	
Acyl. Sample	µl			50

shake for 60 minutes at room temperature

4 x washing

Enzyme Conjugate	µl	100	100	100
------------------	----	-----	-----	-----

shake for 20 minutes at room temperature

4 x washing

Substrate	µl	100	100	100
-----------	----	-----	-----	-----

shake for 15 - 25 minutes at room temperature

Stop Solution	µl	100	100	100
---------------	----	-----	-----	-----

Reading of absorbance at 450 nm