



## Instructions for use

# Serotonin ELISA

Enzyme Immunoassay  
for the Quantitative Determination of Serotonin  
Serum, Plasma and Urine




EA602/96



12 x 8



2 – 8 °C


 DLD Gesellschaft für Diagnostika und medizinische Geräte mbH  
Adlerhorst 15 • 22459 Hamburg • Tel 040-555 87 10 • Fax 040-555 87 111  
Internet: <http://www.dld-diagnostika.de> • E-Mail: [contact@dld-diagnostika.de](mailto:contact@dld-diagnostika.de)



## Contents

1.	Introduction and Principle of the Test.....	5
2.	Precautions.....	6
3.	Storage and Stability.....	6
4.	Contents of the Kit .....	7
5.	Sample Collection.....	9
6.	Preparation of Reagents and Samples .....	9
7.	Test Procedure .....	11
8.	Calculation of the Results .....	12
9.	Assay Characteristics .....	13
10.	Literature .....	15
	Pipetting Scheme.....	16

## Symbols



In Vitro Diagnostic Medical Device



EC Declaration of Conformity



Content



Expiry Date



Batch code



Temperature limitation



Manufacturer



Sufficient for ... determinations



Catalogue number



Consult instructions for use

## Hazard Pictograms



Danger



Warning



## **1. Introduction and Principle of the Test**

Serotonin (5-Hydroxytryptamine), a biogenic amine, is a product of the tryptophan metabolism. It is a well evaluated neurotransmitter of the central nervous system and can be found in high concentrations in the chromaffine cells of the intestinal mucosa, in the platelets and the serotonergic neurones of the brain.

Central-serotonergic neurones influence physiological functions such as sleep and the hormonal and cardio-vascular regulation. Increased serum levels can be found with malignant carcinoid, endogenous depression and schizophrenia.

The assay kit provides materials for the quantitative measurement of derivated serotonin (5-Hydroxytryptamine) in serum, plasma and urine. The derivation is performed during the preparation of the samples. By using the acylation reagent the serotonin is quantitatively derivated into N-acylserotonin.

The competitive Serotonin ELISA kit uses the microtitre plate format. Serotonin is bound to the solid phase of the microtiter plate. Acylated serotonin and solid phase bound serotonin compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. The antibody bound to the solid phase serotonin is detected by anti-rabbit/peroxidase. The substrate TMB / peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase serotonin is inversely proportional to the serotonin concentration of the sample.

## **2. Precautions**

- For in vitro diagnostic use only.
- Disposable gloves and safety glasses should be used.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy but these materials should be handled as potentially infectious.
- Some components of this kit are containing hazardous reagents. These components are marked with the adequate hazard label.

## **3. Storage and Stability**

On arrival, store the kit at 2-8 °C. Once opened the kit is stable until its expiry date. For stability of prepared reagents refer to Preparation of Reagents.

Do not use components beyond the expiration date shown on the kit labels.

Do not mix various lots of any kit component within an individual assay.

## 4. Contents of the Kit

- 4.1. MT-Strips** **STRIPS** 12 strips  
8 wells each, break apart, precoated with serotonin

- 4.2. Standards 1 - 6** **CAL 1-6** 6 vials  
4 ml each, ready for use

Concentrations:

Standard	1	2	3	4	5	6
ng/ml	0	15	50	150	500	2,500

- 4.3. Control 1 & 2** **CON 1 & 2** 2 vials  
4 ml each, ready for use,  
Range: see q.c. certificate

- 4.4. Acylation Buffer** **ACYL-BUFF** 1 vial  
3 ml ready for use, blue coloured



Warning

- 4.5. Acylation Reagent** **ACYL-REAG** 1 vial  
2,5 ml, ready for use, contains DMSO  
and DMF;  
Please note: Acylation Reagent reacts  
with many plastic materials including  
plastic trays. It does not react with  
normal pipette tips and with glass  
devices.




Warning



Danger

- 4.6. Antiserum** **AS** 1 vial  
11 ml, ready for use, colour coded yellow  
Rabbit-anti-N-acetylserotonin

<b>4.7. Enzyme Conjugate</b>	<b>CONJ</b>	1 vial
12 ml, , ready for use Goat anti-rabbit-IgG-peroxidase		 Warning
<b>4.8. Wash Buffer</b>	<b>WASH</b>	1 vial
20 ml, 50x concentrated Dilute contents with dist. water to 1 litre total volume.		
<b>4.9. Substrate</b>	<b>SUB</b>	1 vial
12 ml TMB solution, ready for use		
<b>4.10. Stop Solution</b>	<b>STOP</b>	1 vial
12 ml, ready for use Contains 0.3 M sulphuric acid		
<b>4.11. Reaction plate</b>	<b>ACYL-PLATE</b>	1 piece
for acylation		
<b>4.12. Equalizing Reagent</b>	<b>EQUA-REAG</b>	1 vial
lyophilized dissolve content with 20.5 ml dist. water, dissolve carefully to minimize foam formation		
<b>4.13. Foil</b>	<b>FOIL</b>	3 oieces
Ready for use		

Additional materials and equipment required, but not provided:

- Pipettes (10, 25, 50, 100 and 200 µl) and multipette
- Orbital shaker
- Multichannel pipette or microplate washing device
- Distilled water
- Centrifuge
- Microplate photometer (450 nm)



## 5. Sample Collection

### 5.1. Serum and Plasma

The test can be performed with serum as well as with EDTA plasma. If plasma is to be used care must be taken to get true platelet-free plasma. Otherwise, the Serotonin level has to be related to the number of thrombocytes in the sample. Since the preparation of platelet-free plasma requires special precautions, it is generally recommended to use serum instead of plasma.

Hemolytic and lipemic samples should not be used.

The samples can be stored up to 6 hours at 2 - 8 °C. For a longer storage (up to 6 months) the samples must be frozen at -20 °C

Repeated freezing and thawing should be avoided.

### 5.2. Urine

The total volume of urine excreted during a 24-hours period should be collected and mixed in a single bottle containing 10 - 15 ml of 6 M hydrochloric acid as preservative. Avoid exposure to direct sun light. Determine the total volume and take an aliquot for the measurement. For patients with suspected kidney disorders the creatinine concentration should be tested, too. Urine samples can be stored at -20 °C for at least 6 months.

## 6. Preparation of Reagents and Samples

### 6.1. Microtiter strips **STRIPS**

Before opening the packet of strip wells, allow it to stand at room temperature for at least 10 minutes. After opening, keep any unused wells in the original foil packet with the desiccant provided. Reseal carefully and store at 2-8 °C.

### 6.2. Wash Buffer **WASH**

Dilute the content with dist. water to a total volume of 1,000 ml.

For further use the diluted wash buffer must be stored at 2 - 8 °C for a maximum period of 4 weeks.

### 6.3. Equalizing Reagent **EQUA-REAG**

Dissolve the content with 20.5 ml dist. water, mix shortly and leave on a roll mixer for 30 minutes. Handle carefully in order to minimize foam formation. The reconstituted Equalizing Reagent should be stored frozen at -20 °C and is stable for a minimum of 1 year.

All other reagents are ready for use.

### 6.4. Preparation of Samples (Acylation)

The wells of the reaction plate for the acylation should be used only once. Please mark the respective wells before using.

1. Pipette 20 µl of standard 1 – 6 **CAL 1 - 6**, 20 µl of control 1 & 2 **CON 1 & 2** and 20 µl of serum, 20 µl of urine or 40 µl of plasma into the respective wells of the reaction plate **ACYL-PLATE**. Duplicates are recommended.
2. Pipette 20 µl Acylation Buffer **ACYL-BUFF** into each well.
3. Pipette 200 µl Equalizing Reagent **EQUA-REAG** into each well.
4. Shake plate on an orbital shaker for 10 seconds.
5. Pipette 20 µl Acylation Reagent **ACYL-REAG** into each well, mix immediately.

Please note, that Acylation Reagent reacts with many plastic materials including plastic trays. It does not react with normal pipette tips and with glass devices. Use an Eppendorf multipipette or similar, fill the syringe directly from the vial and add well by well.

6. Incubate for 15 minutes at room temperature (approx. 20 °C) on an orbital shaker. Colour changes to green.
7. Take 20 µl each for the ELISA.

## 7. Test Procedure

### 7.1. Sample Incubation

Pipette 20 µl of acylated Standards 1 to 6, 20 µl of acylated controls and 20 µl of acylated samples into the respective wells of the coated microtiter strips **STRIPS**.

Pipette 100 µl Antiserum **AS** into each well.

Seal plate with foil **FOIL** and incubate for 30 minutes at room temperature (20 – 25 °C) on an orbital shaker.

### 7.2. Washing

Discard or aspirate the contents of the wells and wash thoroughly with 250 µl Wash Buffer **WASH** per well. Repeat the washing procedure 4 times. Remove residual liquid by tapping the inverted plate on clean absorbent paper.

### 7.3. Conjugate Incubation

Pipette 100 µl enzyme conjugate **CONJ** into each well.

Seal plate with foil **FOIL** and incubate for 15 minutes at room temperature on an orbital shaker.

### 7.4. Washing

Repeat step 7.2.

### 7.5. Substrate Incubation

Pipette 100 µl Substrate **SUB** into each well.

Shake on an orbital shaker for 10 seconds and then incubate for 15 ± 5 minutes at room temperature, without shaking, on the table, cover plate with a large box.

### 7.6. Stopping

Pipette 100 µl Stop Solution **STOP** into each well, in the same order and timing as performed in step 7.5.

### 7.7. Reading

Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer.


## 8. Calculation of the Results

On a semilogarithmic graph paper the concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Alternatively, the optical density of each standard and sample can be related to the optical density of the zero standard, expressed as the ratio  $OD/OD_{max}$ , and then plotted on the y-axis.

The concentration of the controls, urine and serum samples can be read directly from this standard curve by using their average optical density.

The read-off values for plasma samples have to be divided by a factor of 1.8.

Typical standard curve:



$y = (A - D) / (1 + (x/C)^B) + D$ :      A      B      C      D      R<sup>2</sup>  
 ○ Std (Standards: Concentration vs MeanValue)      2,205      0,947      119,181      0,159      0,999

Standard	1	2	3	4	5	6
ng/ml	0	15	50	150	500	2,500
nmol/ l	0	85.1	284	851	2,838	14,188

Conversion: Serotonin: 1ng/ ml = 5.675 nmol/l

## 9. Assay Characteristics

### 9.1. Normal Range

The reference ranges given below should only be taken as a guideline. It is recommended that each laboratory establishes its own normal values.

Matrix	Reference range
Serum (Female)	80 - 450 ng/ml
Serum (Male)	40 - 400 ng/ml
Urine	50 - 250 µg/day
Plasma (patelet free)	< 10 ng/ml

### 9.2. Sensitivity

Matrix	Lower Detection Limit	Calculation
Serum, Urine	4.7	OD <sub>Cal1</sub> - 2xSD
Plasma	2.6	OD <sub>Cal1</sub> - 2xSD

### 9.3. Specificity (Cross Reactivity)

Structural related components were tested for possible interference with the antisera against Serotonin used in the ELISA method. The tested compounds were Tryptamine, Melatonin, 5-HIAA, 5-Hydroxy-L-Tryptophan, 5-Methoxytryptamine and L-Tryptophan.

Substance	ED-50-Value (ng/ml)	Cross Reactivity (%)
Serotonin	133	100
Tryptamine	8,700	1.5
5-Methoxytryptamine	56,900	0.23
Melatonin	> 1,000,000	< 0.0133
5-Hydroxy-L-Tryptophan	> 1,000.000	< 0.0133
5-HIAA	> 10,000,000	< 0.00133
L-Tryptophan	> 10,000,000	< 0.00133

### 9.4. Recovery

Increasing amounts of Serotonin were added to a serum, plasma and an urine sample. Each spiked sample was assayed. The analytical recovery

of Serotonin was estimated at different concentrations by using the theoretically expected and the actually measured values.

<b>Matrix</b>	<b>Range (ng/ml)</b>	<b>Mean (%)</b>	<b>Range (%)</b>
<b>Serum</b>	70 - 824	95	85 - 105
<b>Urine</b>	27 - 1085	105	83 - 120
<b>Plasma</b>	62 - 293	96	87 - 102

### 9.5. Linearity

The linearity of the ELISA method was investigated using different dilutions of a serum, plasma and an urine sample. Samples were diluted with distilled water.

<b>Matrix</b>	<b>Range (ng/ml)</b>	<b>max. Dilution</b>	<b>Mean (%)</b>	<b>Range (%)</b>
<b>Serum</b>	60 – 1,203	1 : 20	91	83 - 97
<b>Urine</b>	66 – 1,316	1 : 20	100	96 - 104
<b>Plasma</b>	79 - 395	1 : 5	94	90 - 97

### 9.6. Reproducibility

The reproducibility of the ELISA method was investigated by measuring the intra-assay-coefficients of variation (cv) by repeated measurements of different serum, plasma and urine samples with different Serotonin concentrations.

<b>Matrix</b>	<b>Range (ng/ml)</b>	<b>Intra-Assay Variation cv (%)</b>
Serum	148 – 497	7.3 – 6.9
Urine	93 – 209	6.7 – 6.1
Plasma	163	7.7

## 10. Literature

- Kema, P.; de Vries, E.; Muskiet, F. (2000):  
**Clinical chemistry of serotonin and metabolites**  
Journal of Chromatography B, 747 33–48
- Lechin, F.; van der Dijs, B.; Lechin, A. (2005):  
**Circulating Serotonin, Catecholamines, and Central Nervous System Circuitry Related to Some Cardiorespiratory, Vascular, and Hematological Disorders**  
The Journal of Applied Research Vol. 5, No. 4
- Spaeth, M. (2006):  
**Fibromyalgia Syndrome: The Role of Neurochemicals**  
Primary Psychiatry. 13(9):72-75
- Spivak, B.; Vered, Y.; Graff, E.; et al. (1999):  
**Low Platelet-Poor Plasma Concentrations of Serotonin in Patients with Combat-Related Posttraumatic Stress Disorder**  
BIOL PSYCHIATRY. 45:840–845
- Ernberg, M.; Lundeberg, T.; Kopp, S. (2000):  
**Pain and allodynia/hyperalgesia induced by intramuscular injection of serotonin in patients with fibromyalgia and healthy individuals**  
Pain 85 31±39
- Alvarez, J.; Gluck, N.; Fallet, A.; et al. (1999):  
**Plasma serotonin level after 1 day of fluoxetine treatment: a biological predictor for antidepressant response?**  
Psychopharmacology 143 : 97.101
- Mück-Seler, D.; Pivac, N.; Jakovljevic, M.; et al. (1999):  
**Platelet Serotonin, Plasma Cortisol, and Dexamethasone Suppression Test in Schizophrenic Patients**  
BIOL PSYCHIATRY 45:1433–1439
- Vikenes, K.; Farstad, M.; Nordrehaug, J. (1999):  
**Serotonin Is Associated with Coronary Artery Disease and Cardiac Events**  
*Circulation* August 3, 1999; 483-489
- Dayan, P.; Huys, Q.. (2008):  
**Serotonin, Inhibition, and Negative Mood**  
PLoS Computational Biology February 2008 | Volume 4 | Issue 2 | e4
- Leboyer, M.; Philippe, A.; Bouvard, M.; et al. (1999):  
**Whole Blood Serotonin and Plasma Beta-Endorphin in Autistic Probands and Their First-Degree Relatives**  
BIOL PSYCHIATRY 1999;45:158–163

## Pipetting Scheme

### Sample Preparation

		Standards	Control	Serum, Urine	Plasma
<b>Standard 1- 6</b>	μl	20			
<b>Control 1 &amp; 2</b>	μl		20		
<b>Serum, Urine</b>	μl			20	
<b>EDTA-Plasma</b>	μl				40
<b>Acyl. Buffer</b>	μl	20	20	20	20
<b>Equalizing Reag.</b>	μl	200	200	200	200

Shake plate for 10 seconds on an orbital shaker

<b>Acyl. Reag.</b>	μl	20	20	20	20
--------------------	----	----	----	----	----

Mix immediately and incubate 15 minutes at room temperature on an orbital shaker  
Take 20 μl each for the ELISA

### Pipetting Scheme ELISA

		Standards	Control	Serum, Urine	Plasma
<b>acyl. Standard 1-6</b>	μl	20			
<b>acyl. Control 1 &amp; 2</b>	μl		20		
<b>acyl. Sample</b>	μl			20	20
<b>Antiserum</b>	μl	100	100	100	100

Seal plate with foil

Incubate for 30 minutes at room temperature on an orbital shaker

Wash 4 x

<b>Conjugate</b>	μl	100	100	100	100
------------------	----	-----	-----	-----	-----

Seal plate with foil

Incubate for 15 minutes at room temperature on an orbital shaker

Wash 4 x

<b>Substrate</b>	μl	100	100	100	100
------------------	----	-----	-----	-----	-----

Shake plate for 10 seconds

Incubate for 15 ± 5 minutes at room temperature, without shaking. Cover plate with large box.

<b>Stop Solution</b>	μl	100	100	100	100
----------------------	----	-----	-----	-----	-----

Read absorbance at 450 nm