



Instructions for Use

Serotonin High Sensitive ELISA

(high sensitivity and small sample volume)


Highly Sensitive Enzyme Immunoassay for the
Quantitative Determination of
Serotonin

RUO

For Research Use Only
Not for Use in Diagnostic Procedures

REF EA630/96

 12 x 8

 2 – 8 °C











 DLD Gesellschaft für Diagnostika und medizinische Geräte mbH
Adlerhorst 15 • 22459 Hamburg • Telefon: 040/ 555 87 10 • Fax: 040/ 555 87 111
Internet: <http://www.dld-diagnostika.de> • E-Mail: contact@dld-diagnostika.de

Table of Contents

1	Introduction and Principle of the Test	4
2	Precautions	4
3	Storage and Stability	4
4	Contents of the Kit	5
5	Sample Collection.....	7
6	Preparation of Reagents and Samples.....	8
7	Test Procedure ELISA.....	11
8	Calculation of the Results.....	12
9	Assay Characteristics.....	13
10	Changes to declare.....	14
	Pipetting Scheme Sample Preparation	15
	Pipetting Scheme ELISA.....	16

Symbols

	For Research Use Only		
	Content		Expiry Date
	Lot Number		Store at
	Manufactured by		Sufficient for n determinations
	Catalogue Number		Consult Instructions for Use

Hazard Pictograms

	Danger		Warning
---	--------	---	---------

1 Introduction and Principle of the Test

The Serotonin High Sensitive ELISA provides materials for the quantitative measurement of derivated serotonin in low concentrated samples and for small sample volumes. The derivation is performed during the preparation of the samples. By using the acylation reagent the serotonin is quantitatively derivated into N-acylserotonin.

The competitive ELISA kit uses the microtiter plate format. Derivated serotonin compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. The antibody bound to the solid phase serotonin is detected by anti-rabbit/peroxidase. The substrate TMB/peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase serotonin is inversely proportional to the serotonin concentration of the sample.

2 Precautions

- For research use only. Not for use in diagnostic procedures.
- Disposable gloves and safety glasses should be used.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy but these materials should be handled as potentially infectious.
- Some components of this kit are containing hazardous reagents. These components are marked with the adequate hazard label.

3 Storage and Stability


On arrival, store the kit at 2-8 °C. Once opened the kit is stable until its expiry date. For stability of prepared reagents refer to Preparation of Reagents.



Do not use components beyond the expiration date shown on the kit labels.

Do not mix various lots of any kit component within an individual assay.

All reagents must equilibrate to room temperature before use and be refrigerated immediately after use.

4 Contents of the Kit

MT-Strips	STRIPS	12 strips
8 wells each, break apart, precoated with Serotonin		
Standard	CAL	1 vial
4 ml, concentrated, concentration: 500 ng/ml Dilute concentrate to working concentrations (see 6.1.2.)		
Control 1 & 2	CON 1 & CON 2	2 vials
Each 4 ml, concentrated, Dilute 1:500 (see 6.1.3.), Range: see QC certificate		
Acylation-Reagent	ACYL-REAG	4 vials
Lyoph., dissolve content in 2.5 ml SOLVENT (see 6.1.5.)		
Acylation Buffer	ACYL-BUFF	1 vial
Lyoph., dissolve content with 4 ml distilled water, add 200 µl of Acylation Buffer Concentrate ACYL-BUFF-CONC Mix carefully in order to minimize formation of foam (see 6.1.4.)		
Acylation Buffer Concentrate	ACYL-BUFF-CONC	1 vial
1 ml, concentrated, colour coded yellow		 Warning
Deactivator	DEAC	1 vial
3 ml, ready for use, colour coded blue		
Enzyme Conjugate	CONJ	1 vial
12 ml, ready for use anti-rabbit-IgG-peroxidase		
Wash Buffer	WASH	1 vial
20 ml, concentrated Dilute content with distilled water to 500 ml total volume (see 6.1.6.)		

Substrate 12 ml TMB solution, ready for use	SUB	1 vial
Stop Solution 12 ml, ready for use Contains 0.3 M sulphuric acid	STOP	1 vial
Solvent 6 ml, ready for use, Solvent to dissolve the Acylation Reagent, contains Acetone	SOLVENT	2 vials
		Warning
		Danger
Ascorbic Acid 2 ml, ready for use contains 10% ascorbic acid	ASC-ACID 10%	1 vial
Standard Buffer 50 ml, contains 10 mM PBS (0,9% NaCl), stabilized Before use enrich to 0.1% ascorbic acid (see 6.1.1.)	STD-BUFF	1 vial
Reaction plate For acylation, ready for use	ACYL-PLATE	1 piece
Adhesive foil Ready for use	FOIL	2 pieces

Additional materials and equipment required but not provided:

- Pipettes (10, 20, 25, 50, 100 and 200 µl)
- Orbital shaker
- Multichannel pipette or Microplate washing device
- Microplate photometer (450 nm and 570 – 650 nm)
- Distilled water

5 Sample Collection

The test is intended for small sample volumes, for low concentrated samples (e.g. tissue homogenates, dialysates) and in general for diluted samples.

To protect serotonin against oxidative degradation the samples must contain 0.1% ascorbic acid.

The samples can be stored up to 6 hours at 2 – 8 °C. For a longer storage the samples must be frozen at -20 °C. Repeated freezing and thawing should be avoided.

Different dilution buffers are suitable, but have to be tested beforehand. Evaluation was done with Ringer buffer and PBS (0.9% NaCl). Alternatively, **STD-BUFF** included in the kit can be used. All applied buffers must contain 0.1% ascorbic acid.

For small sample volumes (< 20 µl) a volume correction is necessary. Add your dilution buffer (alternatively prepared Standard Buffer) to correct for volume.

For example:

Sample volume	Volume dilution buffer
1 µl	19 µl
2 µl	18 µl
5 µl	15 µl
10 µl	10 µl
15 µl	5 µl
20 µl	/

6 Preparation of Reagents and Samples

6.1 Preparation of Reagents

6.1.1 Standard Buffer

The Standard Buffer [STD-BUFF] has to be enriched to 0.1 % ascorbic acid prior use: e.g. 50 ml [STD-BUFF] + 0.5 ml [ASC-ACID 10%].

The prepared Standard Buffer should be stored at -20 °C and is stable for a minimum of 1 year.

6.1.2 Standard

The concentration of [CAL] is 500 ng/ml (= 10.000 pg/sample) serotonin.

Dilute [CAL] to obtain working concentrations as follows:

Std 6	100 pg/sample	990 µl Dilution buffer	+	10 µl CAL
Std 5	20 pg/sample	800 µl Dilution buffer	+	200 µl Std 6
Std 4	6.7 pg/sample	933 µl Dilution buffer	+	67 µl Std 6
Std 3	2 pg/sample	980 µl Dilution buffer	+	20 µl Std 6
Std 2	0.67 pg/sample	993 µl Dilution buffer	+	6.7 µl Std 6
Std 1	0 pg/sample	1000 µl Dilution buffer		

Use the same dilution buffer for dilution of standards as present in the samples or used for diluting the samples.

Alternatively, the prepared Standard Buffer can be used for dilution of standards and samples.

All applied buffers (also Standard Buffer) must contain 0.1% ascorbic acid.

Dilution should be done in polypropylene (PP) tubes or polypropylene (PP) microtubes.

The diluted standards must be prepared immediately before use. After use, discard the standards.

6.1.3 Control 1 & 2

The controls **CON 1** & **CON 2** must be diluted 1:500 prior to use:

dil. Control 1 (1:500):	5,000 µl dilution buffer	+	10 µl CON 1
dil. Control 2 (1:500):	5,000 µl dilution buffer	+	10 µl CON 2

Use the same dilution buffer for controls, samples and standards.

6.1.4 Acylation Buffer

Dissolve the content of **ACYL-BUFF** with 4 ml distilled water. Add 200 µl of Acylation Buffer Concentrate **ACYL-BUFF-CONC**.

Mix shortly and leave on a roll mixer for 30 minutes. Handle carefully in order to minimize foam formation. The reconstituted Acylation Buffer should be stored frozen at -20 °C and is stable until the expiration date printed on the vial label.

6.1.5 Acylation Reagent

Dissolve the content of one bottle **ACYL-REAG** with 2.5 ml **SOLVENT** and shake for 5 minutes on an orbital shaker. The Acylation Reagent must be prepared immediately before use. After use, discard the reagent. The kit contains 4 vials allowing a maximum of 4 assay runs.

Please note that the solvent reacts with many plastic materials including plastic trays; solvent does not react with normal pipette tips and with glass devices.

Solvent is volatile and the dissolved Acylation Reagent evaporates quickly. Therefore, please do not use a tray with a big surface together with a multichannel pipette for pipetting the Acylation Reagent. Rather, use an Eppendorf multipipette (or similar device), fill the syringe directly from the vial with dissolved Acylation Reagent and add well by well.

6.1.6 Wash Buffer

Dilute the content of **WASH** with distilled water to a total volume of 500 ml.

For further use, the diluted wash buffer must be stored at 2 - 8 °C for a maximum period of 4 weeks.

All other reagents are ready for use.

6.2 Preparation of Samples (Acylation)

Allow reagents and samples to equilibrate to room temperature.

Determinations in duplicates are recommended.

The wells of the **ACYL-PLATE** can be used only once. Therefore, please mark the respective wells before use to avoid repeated use.

1. Pipette **25 µl prepared Acylation Buffer** into the respective wells of **ACYL-PLATE**.
2. Pipette each **20 µl diluted Standard 1 - 6, diluted Control 1 & 2 and Sample** into the respective wells.
3. Mix the reaction plate for 10 seconds.
4. Pipette **10 µl prepared Acylation Reagent** into each well (colour changes to red) and continue with step 5., **immediately.**

Please note that solvent reacts with many plastic materials including plastic trays; solvent does not react with normal pipette tips and with glass devices.

Solvent is volatile and the dissolved Acylation Reagent evaporates quickly. Therefore, please do not use a tray with big surface together with a multichannel pipette for pipetting Acylation Reagent. Rather, use an Eppendorf multipipette (or similar device), fill the syringe directly from the vial with dissolved Acylation Reagent and well by well.

5. Incubate for 60 minutes at room temperature on an orbital shaker. Avoid direct sunlight.

Do not cover the wells or the plate; leave the plate open on the shaker.

6. Pipette **25 µl DEAC** into each well.
7. Cover the plate with **FOIL**.
8. Incubate for 3 hours at room temperature on an orbital shaker. Avoid direct sunlight.

7 Test Procedure ELISA


1. Pipette **50 µl acylated Standards 1 to 6, Controls and Samples** into the respective wells of the **STRIPS**.
2. Cover the plate with **FOIL** and incubate for 15 – 20 hours (overnight) at 2 - 8 °C.
3. Discard or aspirate the contents of the wells and wash thoroughly with each **300 µl prepared Wash Buffer**. Repeat the washing procedure 3 to 4 times. Remove residual liquid by tapping the inverted plate on clean absorbent paper.
4. Pipette **100 µl CONJ** into each well.
5. Incubate for 60 minutes at room temperature on an orbital shaker.
6. Wash: Repeat step 3.
7. Pipette **100 µl SUB** into each well.
8. Incubate for 20 to 30 minutes at room temperature (20 – 25 °C) on an orbital shaker.
9. Pipette **100 µl STOP** into each well.
10. Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer within 10 minutes.

8 Calculation of the Results

The concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Alternatively, the optical density of each standard and sample can be related to the optical density of the zero standard, expressed as the ratio OD/OD_{max} , and then plotted on the y-axis.

The concentration of the controls and samples (pg/sample) can be read directly from this standard curve by using their average optical density.

Typical standard curve:



$y = ((A - D)/(1 + (x/C)^B)) + D$:

A	B	C	D	R²
2,612	1,089	3,791	0,154	1

○ Std (Standards: Concentration vs MeanValue)

9 Assay Characteristics

9.1 Sensitivity

The lower limit of detection was determined by taking the 2fold standard deviation of the absorbance of the Zero Reference and reading the corresponding value from the standard curve.

Sensitivity : 0.39 pg/sample

9.2 Specificity (Cross Reactivity)

Structural related components were tested for possible interference with the antisera against serotonin used in the ELISA method.

Substance	ED-50-Value (ng/ml)	Cross Reactivity (%)
Serotonin	4.3	100
Tryptamine	1,996	0.22
5-Methoxytryptamine	17,083	0.025
5-Hydroxytryptophan	207,551	0.0021
Melatonin	677,434	< 0.001
5-HIAA	> 2,000,000	< 0.001
L-Tryptophan	> 20,000,000	< 0.0001

9.3 Reproducibility

The reproducibility of the ELISA method was investigated by determination of the intra-assay-coefficients of variation (cv) by repeated measurements of two samples with different serotonin concentrations.

Concentrations in ng/ml

Intra-Assay

Sample	n	Mean Value	sd	cv (%)
1	40	4.7	0.41	8.7
2	40	11.9	0.79	6.6

10 Changes to declare

Version 7: Changes in section 4 are highlighted in gray. Some wordings were revised to provide greater clarity.

Version 6: IFU has been re-formatted.

No changes have been made to components or execution of protocols.

Pipetting Scheme Sample Preparation

		Standards	Controls	Sample
ACYL-PLATE				
Prepared ACYL-BUFF	μl	25	25	25
Dil. Standard 1 – 6	μl	20		
Dil. CON 1 & CON 2	μl		20	
Sample	μl			20

Shake for 10 seconds

Freshly prepared ACYL-REAG	μl	10	10	10
----------------------------	----	----	----	----

Incubate for 60 minutes at room temperature on an orbital shaker
 Do not cover wells or plate, leave the plate open on the shaker
 Avoid direct sunlight

DEAC	μl	25	25	25
------	----	----	----	----

Cover plate with FOIL
 Incubate for 3 hours at room temperature on an orbital shaker.
 Avoid direct sunlight

Pipetting Scheme ELISA

		Standard	Control	Sample
STRIPS				
Acyl. Standard 1 – 6	µl	50		
Acyl. CON 1 & CON 2	µl		50	
Acyl. Sample	µl			50

Cover plate with FOIL
 Incubate for 15 – 20 hours (overnight) at 2-8°C
 3 - 4 x washing

CONJ	µl	100	100	100
------	----	-----	-----	-----

60 minutes incubation at room temperature on an orbital shaker
 3 - 4 x washing

SUB	µl	100	100	100
-----	----	-----	-----	-----

20 - 30 minutes incubation at room temperature on an orbital shaker

STOP	µl	100	100	100
------	----	-----	-----	-----

Reading of absorbance at 450 nm (ref. 570 – 650 nm)